



## ALARM FOR SAFETY CABINETS SPECIFICATIONS



Parc Technologique Delta Sud \* 09340 Verniolle \* France  
Tel + 33 5 61 69 84 43 \* fax + 33 5 61 67 86 03  
info@seat-ventilation.com \* [www.seat-ventilation.com](http://www.seat-ventilation.com)

# GENERAL DESCRIPTION

## Benefits

- Ease of installation and maintenance
- Attractive design
- Digital display in m/s or FPM
- Applications: safety cabinets or hoods

## Features

- Audible and visual alarm
- 3 Push buttons: Fan ON/OFF, light ON/OFF, Mute
- Alarm for past time

## Options

- Alarm relay:  
An alarm relay can be installed optionnaly.
- Battery back up:  
Red LED alarm is still functional up to 6h when unit loses power.
- Digital display in m/s or FPM

## Colors

Plastic fascia is available in 3 colours:white, black and yellow.



## Packing

- Supplied in cardboard box which includes controller circuit board, 'O' seal ring, power supply 220/12V and mounting box.

## International Standards Compliance

- CEM Electromagnetic Compatibility EN 61326:1997 / A1 : 1998 /A2 : 2001 / A3 : 2003 (*Test RC-05-42060-1*) US FCC Part 15 Class B edition 2005 (*Test report RC-05-42060-2-A*).
- RoHS directive
- Laboratory French standard XPX 15206 January. 2005
- CE
- **Labour regulation Nr R4222-13** (Alarm for safety cabinets containing dangerous products)

## Components

- **Circuit board with sensor included**  
Panel mounted circuit board to be installed horizontally onto cabinet with 2X2,8 mm screws; IP55 protection ensured by 'O' seal ring; Face plate is supplied with chemical resistant plastic sticker (horizontal) with control/push buttons on fascia; Mounting box is included.




- **Power supply**  
230Vac to 12Vdc power transformer directly into alarm circuit board.

## OPERATING PRINCIPLE


When the fume extraction fan is running, it causes negative pressure inside the safety cabinet or hood. While normal operation (closed doors) the pressure is intense, the front air velocity is constant. If the doors are opened, the pressure becomes less intense and air velocity decrease.

If an opening is made in the wall of the cabinet, air will enter in at a velocity determine by the same negative pressure that is drawing air in the front of the safety cabinet. The obtained air velocity level measurement by using the numerical sensor is displayed on the fascia by LED illumination or digital read out. An audible alarm will sound in case of opened doors or exhaust fan failure.

The alarm timer count the safety cabinet uptime. It allows to check the filters and installation situation in accordance with recommendation. The setting is made from 0 to 100h and from 5 to 365 days.



# LED & DIGITAL DISPLAY




## SPECIFICATIONS

	LED version	Digital display version
<b>References</b>	<b>819698</b>	<b>819697</b>
<b>Digital display–Led</b>	1 green LED Air flow Ok 1 flashing red LED for Alarm <b>No digital display</b>	1 green LED Air flow OK 1 flashing red LED for Alarm <b>3 digits display with velocity reading</b>
<b>Units</b>	N/A	Meter per second (m/s) or FPM
<b>Display range</b>	N/A	0 - 2.00 m/s
<b>Alarm past time</b>	1 flashing Led 5% remaining time 2 constant Leds past Time	1 flashing Led 5% remaining time 2 constant Leds past Time
<b>Alarm delay</b>	Selectable 15s or 30s	Selectable 15s or 30s
<b>Alarm Indication</b>	1 flashing red LED and audible Buzzer	1 flashing red LED and and audible Buzzer
<b>Alarm mute button</b>		
<b>Light On/Off button</b>		
<b>Fan On/Off button</b>		
<b>Alarm relay</b>	Yes (optional)	Yes (optional)
<b>Battery back up</b>	Yes (optional)	Yes (optional)
<b>Alarm past time</b>	Audible and visual	Audible and visual
<b>Mounting</b>	Mounting box	Mounting box
<b>Calibration</b>	Factory pre-calibaretd @ 0.5m/s . Recalibration possible	Factory pre-calibrated @ 0.5m/s . Recalibration possible
<b>Power supply</b>	12Vdc (power supply included 220v:12v)	12Vdc (power supply included 220v:12v)
<b>Orientation</b>	Horizontal	Horizontal
<b>Dimensions</b>	Front fascia: 210L x 75l x 10H mm Mounting box: 205L x 75l x 25H mm	Front fascia: 210L x 75l x 10H mm Mounting box: 205L x 75l x 25H mm

# DIMENSIONS

## Alarm+mounting box



# OVERALL VIEW

